

Air-Compressed Leg Massager

Relax Your Legs






Air-Compressed Massage helps to...


- Relieves Muscle Soreness
- Relieve Fatigue and Pain
- Improve Blood Circulation
- Relieve Calf Cramps
- Ease Swelling

Operation



- Long press for two seconds to turn on/off
- Blue light: Soft compression with heat
- Red light: Soft compression (no heat)
- Green light: Strong compression with heat
- Self-timing: 15 minutes

Velcro Design—Suitable for any body type



Weight: 390g

Bonus



Extension Velcro

Details

- Input: USB 5V
- Charging port: Micro
- Battery: Li-ion Battery 750 mAh
- Rated Power: 3W
- Charging time: 90 minutes
- Using time: 50 minutes